



Collection jars

IDENTIFICATION

- ▶ Collection jar with snap fitting for coupling to the vacuum regulator PF-38G0-1000 ; PF-38G0-2000; PF-38G0-3000
- ▶ Made of autoclavable transparent polycarbonate
- ▶ Capacity indication by means of silk screened printed graduated scale
- ▶ Cap with outlet hose connector
- ▶ It can be connected to the vacuum regulator using just one hand
- ▶ Conceived to receive themicrobiological filter (cod. KT-2800-0000) which is in the cap to protect the operator
- ▶ Overflow device integrated in the cap
- ▶ Adapter for connection of the jar to G1/2" threaded systems

CE MARKING

Risk class I CE marking as prescribed by Italian legislative decree n. 46 of 24 February "Implementation of Directive 93/42/EEC concerning medical devices" and further modifications

REFERENCE STANDARDS

ISO 10079-3: "Medical suction equipment – Suction equipment powered from a vacuum or pressure source"

HTM 02-01: "Medical gas pipeline systems"


*Validity of the standards is referred to the current year





Collection jars

DIMENSIONS



KEY:

- 1. 500 cc safety jar (cod. PF-4100-1000)
- 2. 150 cc safety jar (cod. PF-4100-2000)
- 3. Microbiological filter (cod. KT-2800-0000 10 pcs box)
- 4. [optional] G 1/2" adapter for jar (cod. PF-380A-0000)

CODE	DESCRIPTION	PACKAGING		
		type	kg	pcs
PF-4100-1000	500 cc Safety jar	bag	0.19	1
PF-4100-2000	150 cc Safety jar	bag	0.14	1
PF-380A-0000	G1/2" Adapter for safety jar	bag	0.03	1
AQ-320V-0713	7 x 13 mm Silicone hose	30 m coil	3.20	-

