



Fastening systems for accessories

ACCESSORIES MOUNTING RAIL

- ▶ Accessories mounting rail in satin finish AISI 304 stainless steel 25 x 10 mm
- ▶ Rail lengths produced to client specifications
- ▶ Finishing with plastic caps at the rail ends
- ▶ Anodized aluminium brackets for fixing the rail to the wall
- ▶ Maximum recommended distance between brackets for wall mounting 80 cm




REFERENCE STANDARDS

ISO 19054: "Rail systems for supporting medical equipment"

*Validity of the standards is referred to the current year

CE MARKING

Risk class I CE marking as prescribed by Italian legislative decree n. 46 of 24 February "Implementation of Directive 93/42/EEC concerning medical devices" and further modifications



KEY:

- 1. Fastening wall rail support
- 2. Stainless steel rail
- 3. Universal rail clamp

FASTENING SYSTEMS FOR ACCESSORIES

CODICE	DESCRIZIONE	IMBALLO	
		type	kg
PF-2900-2001	25 x 10 Stainless steel rail	bag	0.84/m
PF-2900-4000	Wall rail support	bag	0.10
APPLICATIONS FOR PLASTERBOARD WALLS			
AQ-1400-0000	Fixation dowel for plasterboard walls	bag	0.10
MP-14UF-0890	Stainless steel screw TCEI M8x90	bag	0.10
APPLICATIONS FOR CONCRETE WALLS			
AQ-1400-0001	Fixation dowel for concrete walls	bag	0.10
MP-14UF-0850	Stainless steel screw TCEI M 8x50	bag	0.10





Fastening systems for accessories

UNIVERSAL CLAMP FOR ACCESSORIES HOLDER RAIL

- ▶ Made of technopolymer material
- ▶ Suitable for attachment to rail sizes 25 ÷ 35 x 10 mm
- ▶ Facility for mounting 30 x 5 mm or 25 x 5 mm accessory bracket by means of pull-out support
- ▶ Lateral fastening system for flexible hose



WALL SUPPORT SLIDE

- ▶ Made of anodized aluminium
- ▶ Wall fixing by means of TSPEI M5 hex socket screws



REFERENCE STANDARDS

ISO 5359: "Low-pressure hose assemblies for use with medical gases"

ISO 19054: "Rail systems for supporting medical equipment"

HTM 02-01: "Medical gas pipeline systems"

*Validity of the standards is referred to the current year

FASTENING SYSTEM FOR ACCESSORIES

CODE	DESCRIPTION	PACKAGING	
		type	kg
PF-2900-6000	Universal clamp 25-35 mm	bag	0.09
PF-2900-5000	Wall bracket	bag	0.09

