

Gas name rings for terminal units

IDENTIFICATION

- ▶ Front ring for medical gas distribution terminal unit type AFNOR NF S 90-116 and NIST (for Instrument Air)
- ▶ Made of plastic with printed wording
- ▶ Indication of gas name in Italian and English, reference symbol, and CE marking of the manufacturer of the Medical Device according to conventions in compliance with ISO 9170-1
- ▶ Customization on request with client's logo
- ▶ Complete with fastening system for installation of covers PL-2301-1000 and PL-2302-xxxx

REFERENCE STANDARD

ISO 9170-1: "Terminal units for medical gas pipeline systems - Part1: Terminal units for use with compressed medical gases and vacuum"

* Validity of the standards is referred to the current year





CODE	GAS	SYMBOL	COLOR	OUTLET STANDARD	PACKAGING		
					type	kg	pcs
PL-23A0-100V	Vacuum	Vac	Yellow	AFNOR	bag	0.004	1
PL-23A0-10O2	Oxygen	O ₂	White	AFNOR	bag	0.004	1
PL-23A0-10AC	Medical Air	Air	Black - White	AFNOR	bag	0.004	1
PL-23A0-1N2O	Nitrous Oxide	N ₂ O	Blue	AFNOR	bag	0.004	1
PL-23A0-1CO2	Carbon Dioxide	CO ₂	Grey	AFNOR	bag	0.004	1
PL-23N0-1AC8	Instrument Air	Air-800	Black - White	NIST	bag	0.004	1
PL-23A0-100N	Azoto	N ₂	Black	NIST	bag	0.004	1
PL-23N0-10HE	Elio	He	Brown	NIST	bag	0.004	1






Gas name rings for terminal units

DIMENSIONS



INSTALLATION IN BED HEAD PANELS AND PENDANTS



KEY:

1. Color coded ring
2. AFNOR or NIST (for Instrument Air) Outlet
3. Extruded profile AQ-3400-xxxx (on request)

