


Integrated panel solutions (pression regulators groups and shut-off valves)

CABINETS

- ▶ Flush mounting cabinet for the installation of in-line pressure regulators (models PF-05xx-yyyy or PF-37xx-yyyy with integrated area block) e and ball valves for gas shut-off and monitoring of the status of service of the medical gas network distribution
- ▶ Plexiglass lockable transparent door, with breakable windows, for quick access to each line of gas/vacuum
- ▶ Integrated housing of alarm module
- ▶ Provision for connection to earth bonding hub
- ▶ Equipped with disposable panel for wall embedding
- ▶ Possibility of supplying assembled cabinet
- ▶ Facility to customize the cabinet in accordance with the client's requirements



DIMENSIONS



CODE	# VALVE SEATS	DIMENSIONS						PACKAGING		
		FINITE DIMENSIONS			EMBEDDED DIMENSIONS			type	cm	kg
		A mm	B mm	C mm	A1 mm	B1 mm	C1* mm			
CA-0501-3000	1 gas+vac	530	680	30	430	580	100 ± 5	box	57 x 72 x 18	12.8
CA-0502-3000	2 gas+vac	530	960	30	430	860	100 ± 5	box	57 x 100 x 18	17.5
CA-0503-3000	3 gas+vac	530	1240	30	430	1140	100 ± 5	box	57 x 128 x 18	22.2
CA-0504-3000	4 gas+vac	530	1520	30	430	1420	100 ± 5	box	57 x 156 x 18	26.9

* Recess depth, follow indicated tolerances

