




## Surface mount installation systems for terminal units





### IDENTIFICATION

- ▶ Surface mount wall outlet box for installation of medical gas and AGSS terminal units (model for systems with centralized blower)
- ▶ Equipped with wall box cover with bas-relief moulding for gas identification sticker
- ▶ Construction in white ABS with threaded inserts for fastening terminal units
- ▶ Modular design to allow multiple boxes to be installed in series

### SURFACE MOUNT BOX

CODE	DESCRIPTION	PACKAGING		
		 type	 kg	 pcs
PF-2100-4000	Surface mount box for terminal units	unboxed	0.14	1

### GAS IDENTIFICATION STICKERS

CODE	GAS	SYMBOL	COLOR	PACKAGING	
				 type	 pcs
ET-2100-000V	Vacuum	Vac	Yellow	bag	1
ET-2100-000O2	Oxygen	O <sub>2</sub>	White	bag	1
ET-2100-00AC	Medical Air	Air	Black - White	bag	1
ET-2100-0N2O	Nitrous Dioxide	N <sub>2</sub> O	Blue	bag	1
ET-2100-0CO2	Carbon Dioxide	CO <sub>2</sub>	Grey	bag	1
ET-2100-0AC8	Instruments Air	Air-800	Black - White	bag	1
ET-2100-00EV	Anaesthetic gas scavenging	AGSS	Magenta	bag	1











## Surface mount installation systems for terminal units



### DIMENSIONS



Gas identification sticker ET-2100-0xxx



#### KEY:

1. Surface standard box PL-2100-4000
2. AFNOR / NIST configuration
3. Gas colour coded ring PL-23A0-1xxx

