



Surface mount installation systems for terminal units



IDENTIFICATION

- ▶ Surface mount wall outlet box for installation of medical gas distribution and AGSS terminal units
- ▶ Equipped with wall box cover with bas-relief moulding for gas identification sticker
- ▶ Construction in white ABS with threaded inserts for fastening terminal units
- ▶ Modular design to allow multiple boxes to be installed in series

SURFACE MOUNT BOX

CODE	DESCRIPTION	PACKAGING		
		type	kg	pcs
PF-2100-4000	Surface mount box for terminal units	unboxed	0.14	1

GAS IDENTIFICATION STICKERS





CODE	GAS	SYMBOL	COLOR	PACKAGING	
				type	pcs
ET-2100-000V	Vacuum	Vac	Yellow	bag	1
ET-2100-00O2	Oxygen	O ₂	White	bag	1
ET-2100-00AC	Medical Air MA4	Air	Black - White	bag	1
ET-2100-0N2O	Nitrous Oxide	N ₂ O	Blue	bag	1
ET-2100-0CO2	Carbon Dioxide	CO ₂	Grey	bag	1
ET-2100-0AC8	Air for Driving Surgical Tools MA7	Air-800	Black - White	bag	1
ET-2100-00EV	Anaesthetic gas scavenging	AGSS	Magenta	bag	1







Surface mount installation systems for terminal units

DIMENSIONS



Gas identification sticker ET-2100-0xxx



KEY:

1. Surface standard box PL-2100-4000
2. AFNOR / NIST configuration
3. Gas colour coded ring PL-23A0-1xxx

